

MISSION CONTROL™ Blood Gas and Electrolyte Control

Level 3

PN: DD-92003D

LOT: MC-1009

Exp: 2009/11

Expected Ranges Chart

Blood Gas/ISE Analyzers	pH		pCO ₂ mmHg		pO ₂ mmHg		Na ⁺ mmol/L			K ⁺ mmol/L			Ca ⁺⁺ mmol/L			Cl ⁻ mmol/L			Li ⁺ mmol/L			tCO ₂ mmol/L			
	Mean	Min Max	Mean	Min Max	Mean	Min Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	
AVL Scientific																									
945, 947	7.605	7.575 - 7.635	22	18 - 26	151	139 - 163																			
990, 995	7.605	7.575 - 7.635	22	19 - 26	147	136 - 158																			
Compact Series																									
982, 983, 985	7.615	7.573 - 7.645	22	19 - 25	144	136 - 151																			
986									162	157 - 167	7.1	6.6 - 7.6			119	114 - 124	2.05	1.86 - 2.26							
984, 987	7.625	7.595 - 7.665							162	157 - 167	7.1	6.6 - 7.6	0.65	0.51 - 0.80	119	114 - 124					26	22 - 30			
OMNI	7.615	7.585 - 7.645	22	19 - 26	136	125 - 146	160	155 - 165	7.0	6.5 - 7.5	0.58	0.49 - 0.68	120	115 - 125											
9110, 9140	7.625	7.595 - 7.655							160	155 - 165	7.0	6.5 - 7.5	0.64	0.50 - 0.79											
9120, 9130									157	152 - 162	7.0	6.5 - 7.5			120	115 - 125									
9180, 9181									154	149 - 159	7.0	6.5 - 7.5	0.58	0.49 - 0.68	118	113 - 123	2.19	2.00 - 2.40							
Ciba-Corning/Bayer																									
238	7.64	7.61 - 7.67	22	18 - 26	130	122 - 137																			
248	7.635	7.605 - 7.665	22	18 - 26	131	123 - 138																			
348	7.645	7.615 - 7.675	22	19 - 25	141	134 - 148	163	158 - 168	7.0	6.5 - 7.5	0.53	0.44 - 0.63	117	112 - 122											
278	7.645	7.615 - 7.675	22	19 - 25	141	134 - 148																			
280	7.645	7.615 - 7.675	22	19 - 25	141	134 - 148																			
288	7.645	7.615 - 7.675	21	18 - 25	142	135 - 148	162	157 - 167	7.0	6.5 - 7.5	0.53	0.44 - 0.63	117	112 - 122											
664							162	157 - 167	6.8	6.3 - 7.3			118	113 - 123							26	22 - 30			
614, 644							162	157 - 167	6.8	6.3 - 7.3			118	113 - 123											
634	7.65	7.62 - 7.68									0.56	0.42 - 0.71													
654							162	157 - 167	6.9	6.4 - 7.4					2.05	1.86 - 2.26									
800 Series*	7.665	7.635 - 7.695	22	18 - 26	137	128 - 145	162	157 - 167	6.9	6.4 - 7.4	0.50	0.36 - 0.65	118	113 - 123											
Rapid 400, 405	7.665	7.635 - 7.695	22	18 - 26	137	128 - 145	162	157 - 167	6.9	6.4 - 7.4	0.50	0.36 - 0.65	118	113 - 123											
IL																									
1304, 1306, 1312	7.615	7.585 - 7.645	23	20 - 27	139	131 - 147																			
BG3	7.625	7.595 - 7.655	22	19 - 25	139	131 - 147																			
BGE	7.635	0.605 - 7.665	23	20 - 26	138	131 - 145	159	154 - 164	6.6	6.1 - 7.1	0.59	0.50 - 0.69	114	109 - 119											
1610, 1620	7.635	7.605 - 7.665	21	18 - 25	139	129 - 148																			
1630, 1640, 1650	7.635	7.605 - 7.665	21	18 - 24	137	129 - 145	159	154 - 164	6.6	6.1 - 7.1	0.58	0.49 - 0.68	114	109 - 119											
Synthesis 10, 15, 20, 25	7.615	7.585 - 7.645	23	20 - 27	139	131 - 148	158	153 - 163	6.7	6.2 - 7.2	0.58	0.49 - 0.68	114	109 - 119											
Gem Premier	7.635	7.605 - 7.665	21	18 - 24	142	133 - 151	159	154 - 164	6.6	6.1 - 7.1	0.58	0.49 - 0.68													
GEM 3000	7.635	7.605 - 7.665	21	18 - 24	142	133 - 151	159	154 - 164	6.6	6.1 - 7.1	0.58	0.49 - 0.68													
NOVA																									
Electrolyte Systems	7.652	7.622 - 7.682					165	160 - 170	7.7	7.2 - 8.2	0.55	0.41 - 0.70	117	112 - 122	2.08	1.89 - 2.29	28	24 - 32							
Stat Profile 1-5	7.635	7.605 - 7.665	22	18 - 26	142	137 - 147	163	158 - 168	6.9	6.4 - 7.4	0.58	0.49 - 0.68	117	112 - 122											
Radiometer																									
ABL 3, 30	7.645	7.615 - 7.675	21	18 - 25	146	139 - 153																			
ABL 300, 330	7.645	7.615 - 7.675	22	19 - 26	144	136 - 151																			
ABL 5	7.63	7.60 - 7.66	21	18 - 25	142	133 - 151																			
ABL, 50, 500, 510, 520	7.635	7.605 - 7.665	22	18 - 26	143	134 - 151																			
ABL 505	7.635	7.605 - 7.665	21	18 - 25	143	134 - 151	160	155 - 165	6.8	6.3 - 7.3	0.59	0.49 - 0.69													
ABL 555	7.635	7.605 - 7.665	21	18 - 25	143	134 - 151	160	155 - 165	6.8	6.3 - 7.3	0.59	0.49 - 0.69													
ABL 600, 610, 620	7.635	7.605 - 7.665	21	18 - 26	143	134 - 151	160	155 - 165	6.7	6.2 - 7.2	0.59	0.50 - 0.69	112	107 - 117											
EML-100							159	154 - 165	6.7	6.2 - 7.2	0.60	0.51 - 0.70	112	107 - 117											
ABL 70,77	7.650	7.620 - 7.680	22	18 - 26	130	122 - 138	157	152 - 162	6.7	6.2 - 7.2	0.60	0.50 - 0.70	115	110 - 120											
ABL 700 Series**	7.635	7.605 - 7.665	21	18 - 25	143	134 - 151	160	155 - 165	6.8	6.3 - 7.3	0.59	0.50 - 0.69	112	107 - 117											
i-STAT																									
BG, E+	7.665	7.635 - 7.695	22	18 - 26	137	128 - 145	160	155 - 165	6.8	6.3 - 7.3	0.59	0.49 - 0.69	112	107 - 117											
Medica, iLyte, Menarini																									
EasyLyte _{Na/K, Na/K/Cl, Na/K/Li, Na/K/Cl/Li}							155	150 - 160	6.6	6.1 - 7.1			118	113 - 123	2.19	1.90 - 2.50									

*Includes 840, 845, 850, 855, 860 Analyzers

**Includes 705, 710, 715, 720, 725

Diamond Diagnostics recognizes all trademarks and copyrights referenced herein.

Manufacturer and Product Information

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
For Technical Assistance call:
Diamond Diagnostics Technical Services at 1-508-429-0450

Intended Use: MISSION CONTROL™ Blood Gas and Electrolyte Control is an assayed quality control material intended for monitoring the measurements of pH pCO₂, pO₂ in blood gas analyzers and sodium, potassium, chloride, lithium, ionized calcium and total carbon dioxide in ISE electrolyte analyzers.

Product Description: This control material is provided for monitoring analyzer performance. It is packaged in sealed glass ampules, each containing approximately 2 ml of solution. Ampules are packaged 10 per tray with each box containing 3 trays, for a total of 30 ampules per box.

Active Ingredients: MISSION CONTROL™ is a buffered solution of electrolytes (Na⁺, K⁺, Cl⁻, Ca⁺⁺, Li⁺, HCO₃⁻/CO₃⁻²). It has been equilibrated with specific levels of CO₂, O₂, and N₂. This control contains no human-based materials.

For in vitro diagnostics use.

Directions for Use

The control should be brought to a temperature of 20-23°C before use (see instructions regarding Expected Ranges). Allow at least four (4) hours for ampules to equilibrate to this temperature prior to testing.

For pH/blood gas values, the control should be analyzed within one (1) minute of opening. For electrolyte measurements, this product is stable for up to one (1) hour after opening.

Follow the procedures listed below:

1. Before use, hold the ampule at the top and bottom (with forefinger and thumb) and shake 15-20 times (about 10 seconds) to mix the solution. Tap the ampule to restore the liquid to the bottom on the ampule.
2. Open the ampule by snapping off the tip at the score. Use gauze, tissue, gloves, or an appropriate ampule opener to protect fingers from cuts.
3. Immediately introduce the liquid from the ampule to the analyzer. Follow the manufacturer's instructions for sampling a control material. Depending on the sampling procedure chosen, the following instructions apply:
 - a. Direct Aspiration: Sample the control directly from the ampule.
 - b. Syringe Transfer:
 - i. Use a clean, gas-tight syringe attached to a clean, blunt syringe needle (if available).
 - ii. Prime the syringe by slowly aspirating a small amount (0.2-0.3 ml) of solution from the ampule.
 - iii. Discard this liquid, leaving the dead space of the syringe filled with the control.
 - iv. Aspirate the control from the ampule into the primed syringe. Be careful that air is not drawn in with the liquid. Expel 1 to 2 drops, detach the needle and immediately inject the control into the analyzer sample port.
 - c. Ampule Injector/Dispenser: Assemble and fill the ampule injector following the manufacturer's instructions. Expel one or two drops to rinse the outlet tip and inject the control into the analyzer sample port.
 - d. Capillary Mode:
 - i. Install the appropriate adapter for micro sampling onto the instrument.
 - ii. Sample the contents of the ampule following the recommendations of the instrument manufacturer. Be certain to keep the sampling tip of the adapter below the surface of the liquid during aspiration.

Limitations**Limitation:**

1. This control is sensitive to many instrument related factors that affect analytical results. Because it is not a blood-based material, it may not detect certain malfunctions, which would affect the testing of blood.
2. This product is intended for use as a quality control material and can assist in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.

Storage:

Store at 18-25°C. Avoid freezing and exposure to temperatures greater than 30°C. You may also store at 4-25°C without adverse effect.

Expected Ranges:

The values for each control analyte on the enclosed Expected Ranges Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range for these ampules when tested at 23°C. (Note: pO₂ values will vary inversely by about one percent (1%) per degree C that the temperature of the ampules varies from 23°C.

The Expected Ranges are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own expected values and control limits. The mean value established should fall within the Expected Ranges shown on the chart.

