
Mission Trinity™ B

Blood Gas, Electrolyte, Metabolite, CO-Oximeter (Siemens/Bayer) Control - Level 2

English DEUTSCH FRANÇAIS ESPAÑOL PORTUGUÉS CHINESE Русский

Intended Use:
MISSION Trinity B™ Blood Gas, Electrolyte, Metabolite, CO-Oximeter Control is an assayed quality control material used for monitoring the performance of blood gas, electrolyte, metabolite, and CO-Oximeter (Bayer) instrumentation for the analytes and analyzers listed on the Expected Values Chart.

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Product Description:
This control material is provided in three (3) distinct levels of pH, pCO₂, pO₂, Na+, K+, Cl-, Ca⁺⁺, glucose, lactate, tHb, O₂Hb, COHb, methHb & HHb covering the significant range of the instrument performance. It is packaged in sealed glass ampules, each containing 1.8 ml of solution.

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Active Ingredients:
MISSION Trinity B™ is a buffered solution of electrolytes, glucose, lactate and dyes. It has been equilibrated with specific levels of CO₂, O₂, and N₂. This control contains no human or biological materials.

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Directions for Use:
Immediately introduce the liquid from the ampule to the analyzer, following the instrument manufacturer's instructions for sampling a control material. Use direct aspiration, syringe transfer, or capillary mode techniques.

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Limitation:
1. This control is sensitive to many instrument related factors that affect analytical results. Because it is not a blood-based material, it may not detect certain malfunctions, which would affect the testing of blood.

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2. This product is intended for use as a quality control material and can assist in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.

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Storage:
Store at 2-8 °C. Avoid freezing and exposure to temperatures greater than 30 °C.

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Expected Ranges:
The values for each control analyte on the enclosed Expected Ranges Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range for these ampules when test at 25 °C. (Note: pO₂ values will vary inversely by about one percent (1%) per degree C that the temperature of the ampules varies from 25 °C.

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The Expected Ranges are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own expected values and control limits. The mean value established should fall within the Expected Ranges shown on the chart.

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