## English

## Intended Use:

MISSION COMPLETE™ Linearity Control are assayed materials used for confirming the calibration and linearity of blood gas, electrolyte, and metabolite instruments for the analytes and analyzers listed on the Expected Values Chart.

## Product Description:

This control material is provided in five (5) distinct levels of pH, pCO<sub>2</sub>, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, covering the significant range of the instrument performance. It is packaged in sealed glass ampoules, each containing 1.8 ml of solution. Ampoules are packaged in kits containing four (4) ampoules of each level.

## Active Ingredients:

MISSION COMPLETE™ is a buffered solution of electrolytes, glucose, and lactate. It has been equilibrated with specific levels of CO<sub>2</sub>, O<sub>2</sub>, and N<sub>2</sub>. This control contains no human or biological materials.

## Directions for Use:

Immediately introduce the liquid from the ampule to the analyzer, following the instrument manufacturer's instructions for sampling a control material. Use direct aspiration, syringe transfer, or capillary method techniques.

## Limitation:

1. The Linearity Control is sensitive to many instrument related factors that affect analytical results. Because it is not a blood-based material, it may not detect certain malfunctions that would affect the testing of blood.

2. This product is intended for use in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.

## Storage:

The expiration date stated on the Linearity Control packaging is for product stored at 2-8°C. The product may also be stored at room temperature (up to 25°C) for six (6) months, provided the labeled expiration date is not exceeded. Avoid exposure to freezing and temperatures greater than 30°C.

## Expected Ranges:

The values for each analyte on the enclosed Expected Values Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range for ampoules that are at 25°C when tested. (Note: pO<sub>2</sub> values will vary inversely by about one percent (1%) per degree Celsius that the temperature of the ampoules varies from 25°C.)

The Expected Ranges are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own acceptance criteria.

For *In Vitro* Diagnostic Use  
In Vitro Diagnostikum  
Usage *In Vitro*

Para Uso Diagnóstico *In Vitro*  
Utivar *Aparato* *In Diagnostico* *In Vitro*  
*In Vitro* diagnóstico  
仅供体内诊断使用

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