Product And Manufacturer Information

Contents: 1 - Chloride Electrode
1 – O-ring

Intended Use: Use with Cobas ISE Modules

Storage: Store the Chloride electrode at room temperature (18-25°C) until you are ready to use it in your system.

Precautions: Please check that the electrode has Fill Solution and no air bubbles in the fill chamber near the membrane sample path area prior to installation.

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
For Technical Assistance call:
Diamond Diagnostics Technical Services at 1-508-429-0450

Installation

- 1 Remove the green tape.
- 2 Check that the electrode has Fill Solution and the inner electrode touches or is immersed in the fill solution.
- 3 Check for air bubbles in the fill chamber near sample path. If there are any, tap the electrode body to dislodge air bubbles.
- 4 Verify that the O-ring is in place on the electrode body.
- 5 Install electrode in your system.
- 6 Allow the electrode to equilibrate for at least 15 minutes.
- 7 Perform a calibration.
- 8 Analyze a quality control sample.

Warranty

Warranty is calculated 3 months from the date of shipment from Holliston, MA.

Limitations

Upon installation of this electrode, you may experience a higher Chloride recovery for patient samples than your previous electrode. Diamond Diagnostics Inc. utilizes a more sensitive membrane material than that of Roche that may cause differences in recovery of serum samples only. The difference in recovery may vary by laboratory due to instrument, reagent lot and electrode variability.

Should you determine a difference in Cl recovery, you can utilize the **Conversion Factor** option available inside your Cobas system software. To locate it proceed as listed below where you are able to input the conversion and offset factor you choose.

Program >>>> Test >>>> CI >>>> Calculation

The Roche Cobas Mira Plus Manual, **Data Analysis Section 7, page 7.42** gives the following example:

The result of the control and sample can be linearly converted, using the CONVERSION FACTOR and OFFSET programmed in the test file.

$$C = [CEND \cdot FCORR] + OCORR$$

C Corrected concentration of control or sample
FCORR Conversion factor (test parameter)

CEND Concentration of control or sample
OCORR Conversion offset (test parameter)

As a further example, see the table below:

OEM Cl Value	Diamond Cl Value	Conversion Factor	Offset factor	Corrected Value
100	75	1.333	0	100
100	105	0.95	0	100
100	100	1.00	0	100
100	118	0.847	0	100

Note that the values listed for Diamond are examples only and do not reflect actual results. Values are examples only to show how the conversion factor would work.

Should have any additional questions, please feel free to contact Technical Services at: **1-508-429-0450**.

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The names of various manufacturers, their instruments, and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Diamond Diagnostics expressly disclaims any affiliation with them or sponsorship by them.