Manufacturer and Product Information

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
For Technical Assistance call:
Diamond Diagnostics Technical Services at 1-508-429-0450

Intended Use: Creatine Reagent is for the quantitative determination of creatinine in serum, plasma or urine on the Beckman Synchron CX® & CX® Delta Systems.

Summary And Principle: This product is intended to serve as a functional equivalent to pre-existing material distributed by the Original Equipment Manufacturer (OEM).

Clinical Significance: Creatinine measurements are used in the diagnosis and treatment of renal diseases, in renal dialysis, and as a calculation basis for measuring other urine analytes.

Methodology: Creatinine is determined by mixing a sample with the alkaline picric reagent, creatinine from the sample combines with the reagent to produce a red-colored complex.
Creatinine + Picric Acid → Creatinine-Picrate Complex (red)
Absorbance readings are taken at both 520 nm and 560 nm at 25.6 seconds after sample addition. The differential absorbance is a direct measure of the concentration of creatinine in the sample.

Reagents: Creatine Reagent Kit, BK-443340D

Containing: Alkaline Buffer – 3 x 1600mL is an aqueous solution containing of 0.188 mol/L NaOH, Buffers buffered with sodium borate and sodium phosphate, surfactant, and preservative.
Picric Acid Solution - 3 x 400 mL - containing 0.05 M Picric Acid.

For *in vitro* diagnostic use only

Cautions: Exercise normal laboratory precautions. If contact occurs with skin, rinse affected area with water. If contact with eyes occurs, immediately rinse with copious amount of clean water or eye rinse. In cases of accidental ingestion, contact a physician immediately. Protective clothing and gloves should be worn when handling the picric solution and working solution as the picric acid stains both skin and clothing.

Storage & Stability: The alkaline buffer and picric solution are stable until the expiration date listed on the bottle. The product should not be used beyond this date. The combined working reagent is stable for 30 days from the time of mixing if stored properly at ambient temperatures. Store upright at room temperature, 18-25 °C.
Do not freeze or refrigerate the solution or working reagent. Precipitation may occur at low temperatures. The precipitate will dissolve when warmed at 37°C. Do not use until all of the precipitate has dissolved. Mix by inversion before use.

Procedure

Procedure: Intended to serve as a direct replacement to pre-existing materials distributed by the OEM. Do not prepare working reagent until ready to use. The working solution is prepared prior to installing on the instrument by: pouring a bottle of Picric Acid solution into a bottle of the Alkaline Buffer Reagent, capping and mixing by gentle inversion. For a detailed description of the use of this reagent, refer to the Instrument's Operator Manual.

Quality Control: Diamond Diagnostics suggests the use of commercially available control material with results assayed for the instrument used. Controls should be run at Normal and Abnormal levels. Diamond Diagnostics suggests measuring controls before patient samples are run and following instrument maintenance.

Limitations

Serum, Plasma or Urine are the sample of choice, whole blood samples are not recommended. Compatible anticoagulants or Chemical additives with this method are:

Anticoagulants / Acceptable Level		Anticoagulants / Acceptable Level		Anticoagulants / Acceptable Level	
Sodium Citrate*	3.5 mg/mL	Ammonium Heparin	45 Units/mL	Sodium Heparin	45 Units/mL
EDTA	4.0 mg/mL	Lithium Heparin	45 Units/mL	Lithium Iodoacetate	1.5 mg/mL

*Specimens collected in a tube containing liquid Sodium Citrate will exhibit a decrease in value to dilution effect.
Lipemic samples >+3 should be ultracentrifuged and the analysis performed on the supernatant.

Refer to Literature References 2 and 3 for other interferences caused by drugs and disease.

If the instrument fails calibration or controls do not measure within acceptable range when Diamond Diagnostics products are used, Diamond Diagnostics suggests the following:
Verify that the reagents and internal calibrators used to standardize the instrument are correct for the instrument, have adequate expiration, and do not contain visually evident contamination.
Follow the procedures delineated within the Operator's Manual listed under Troubleshooting.
Ensure that all appropriate Maintenance Procedures, as listed in the Operator's Manual, have been performed.
If problems still exist, contact Diamond Diagnostics' Technical Service Department.

Useable Range:	Creatine was found to have the following useable range: 0.6 to 25 mg/dL for Serum or Plasma 10 to 400 mg/dL in Urine* *Samples with concentrations exceeding the high end of the useable range should be diluted with deionized water or saline and reanalyzed.						
Functional sensitivity:	Evaluated on dilutions of serum samples made from a starting serum of an approximately concentration of 0.9 mg/dL; dilutions were 1:3, 1:5, 1:11 and a 0 mg/dL. Samples were tested as 4 samples per run over 5 calibrated runs. The lowest level where a CV of less than 20% was achieved was with the dilution at an expected value of 0.9 mg/dL BUN which recovered as: 0.6 mg/dL with Diamond reagent 0.8 mg/dL with Beckman reagent						
Reference Ranges^{3,4,5}:	The values given as reference ranges are intended only as a guide rather than as absolute values. No account is taken for age, sex, or dietary differences that may exist. Each lab should establish a reference range based on their patient population. <table border="1"> <thead> <tr> <th>Sample Type</th><th>SI Units</th></tr> </thead> <tbody> <tr> <td>Serum or Plasma</td><td>0.6 to 1.2 mg/dL</td></tr> <tr> <td>Urine (timed)</td><td>0.8 to 2g/24hrs</td></tr> </tbody> </table>	Sample Type	SI Units	Serum or Plasma	0.6 to 1.2 mg/dL	Urine (timed)	0.8 to 2g/24hrs
Sample Type	SI Units						
Serum or Plasma	0.6 to 1.2 mg/dL						
Urine (timed)	0.8 to 2g/24hrs						

Performance Characteristics

Performance data was collected on CX®3 Delta and CX®9ALX Systems

Method Comparison: Diamond Creatinine Reagent was compared to Beckman Reagent using spiked or diluted serum samples run in triplicate and tested with each reagent, Diamond Creatinine Reagent and Beckman Creatinine Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained:
Diamond = $0.988 \times \text{Beckman} + 0.634 \text{ mg/dL}$
 $r^2 = 0.999$ $n=45$ range = 10 to 400 mg/dL $S_{(y,x)} = 3.17 \text{ mg/dL}$

Recovery to Expected Values: Was evaluated for each matrix; serum and urine. Dilutions of the respective matrices were made and measured with Diamond and Beckman reagent. % Recovery = (Measured/expected) x100 was calculated for both Diamond and Beckman. Diamond and Beckman exhibited similar recoveries across the range of values in all matrices. See table 1 below.

Precision: Commercially available control sera and urine controls were assayed using both Diamond and Beckman Creatinine Reagents on Synchron CX® Delta System. Samples were run for 20 days, 2 runs per day, 2 observations per run. See table 2 below.

1				
Matrix	Range of Conc. Expected, mg/dL	Reagent	Range of average % Recovery	Overall Mean Recovery
Serum	9.04 – 0.20 mg/dL	Diamond	85.7 – 105.3	98.5
		Beckman	93.4 – 114.3	104.2
Urine	400 – 10 mg/dL	Diamond	87 – 102.1	96
		Beckman	83 – 102.9	96.2

2				
	N	Mean mg/dL	within run sd	Total sd
Serum Control 1	80	1.6	0.07	0.14
Serum Control 2	80	6.9	0.10	0.62
Urine Control 1	80	89	0.7	10.9
Urine Control 2	80	217	2.2	9.5

Literature References

- Jaffe, M.Z., Physiological Chemistry 10:391 (1886)
- Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
- Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

SYNCHRON CX®3 Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Diamond Diagnostics Technical Services at 1-508-429-0450

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1

Matrix	Range of Conc. Expected, mg/dL	Reagent	Range of average % Recovery	Overall Mean Recovery
Serum	9.04 – 0.20 mg/dL	Diamond	85.7 – 105.3	98.5
		Beckman	93.4 – 114.3	104.2
Urine	400 – 10 mg/dL	Diamond	87 – 102.1	96
		Beckman	83 – 102.9	96.2

2

	N	Mean mg/dL	within run sd	Total sd
Serum Control 1	80	1.6	0.07	0.14
Serum Control 2	80	6.9	0.10	0.62
Urine Control 1	80	89	0.7	10.9
Urine Control 2	80	217	2.2	9.5

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1. Jaffe, M.Z., Physiological Chemistry 10:391 (1886)
2. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
3. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
4. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Serviciul Technic Diamond Diagnostics 1-508-429-0450

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Literatură de referență	
1.	Jaffe, M.Z., Physiological Chemistry 10:391 (1886)
2.	Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2 nd Edition, W.B. Saunders Company, Philadelphia (1994)
3.	Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20 th Edition, W.B. Saunders Company, Philadelphia (2001)
4.	SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA
SYNCHRON CX®3 Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA	

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