Manufacturer and Product Information

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
For Technical Assistance call:
 Diamond Diagnostics Technical Services at 1-508-429-0450


Intended Use: Diamond BUN Reagent, in conjunction with Synchron CX® System Calibrators, is intended for the quantitative determination of urea nitrogen in serum, plasma on Beckman Synchron CX® & CX® Delta Systems and for the determination of urine on CX® Delta Systems.

Summary And Principle: This product is intended to serve as a functional equivalent to pre-existing material distributed by the Original Equipment Manufacturer (OEM).

Clinical Significance: BUN measurements are used in the diagnosis and treatment of certain renal and metabolic diseases. BUN is regulated by the metabolism of proteins and by the renal excretion of Urea.

Methodology: Urea nitrogen concentration is determined by an enzymatic conductivity rate method employing a Beckman Conductivity Electrode. Electronic circuits determine the rate of increase in conductivity, which is directly proportional to the concentration of BUN in the sample^{5,6}.

Chemical Reaction Scheme: When a 10 micro liters sample is injected into a reaction cup containing the urease reagent, the urea in the serum undergoes changes according to the following reaction:

The result to the reaction equation is the conversion of the non-ionic species (urea) to one, which is ionic (ammonium carbonate). During reaction, the timed rate of increase of solution conductivity is directly proportional to the concentration of urea then present in the reaction cup.

Reagents: Bun Reagent, BK-443350D
Containing: Bun Reagent – 500mL is an aqueous solution containing of 120 U/mL Jack Bean Urease, salts, buffer and preservative. Also required on the CX® for the determination of BUN: 1-10 liter bottle of Wash Solution (BK-443335D –reconstituted for 250mL concentrate); Synchron Calibration Standards; Appropriate Controls.

For in vitro diagnostic use only

Cautions: Exercise normal laboratory precautions. If contact occurs with skin, rinse affected area with water. If contact with eyes occurs, immediately rinse with copious amount of clean water or eye rinse. In cases of accidental ingestion, contact a physician immediately.

Storage and Stability: BUN Reagent stored at 2°C to 8°C is stable until the expiration date indicated on each bottle. BUN Reagent stored at room temperature (ie as when on the instrument) is stable for 30 days until expiration date, if sooner. Do not use if frozen.

Procedure

Reagent Preparation: Prior to use, allow the reagent to warm to ambient temperature. An 18 – 25°C water bath may be used to warm reagent.

Quality Control: Diamond Diagnostics suggests the use of commercially available control material with results assayed for the instrument used. Controls should be run at Normal and Abnormal levels. Diamond Diagnostics suggests measuring controls before patient samples are run and following instrument maintenance.

Limitations

Serum, Plasma or Urine are the sample of choice, whole blood samples are not recommended^{5,6}. Compatible anticoagulants or Chemical additives with this method ^{5,6}:

Anticoagulants / Acceptable Level	Anticoagulants / Acceptable Level	Anticoagulants / Acceptable Level
Sodium Citrate* 3.5 mg/mL	Ammonium Heparin 45 Units/mL	Sodium Heparin 45 Units/mL
EDTA 4.0 mg/mL	Lithium Heparin 45 Units/mL	Lithium Iodoacetate 1.5 mg/mL

*Specimens collected in a tube containing liquid Sodium Citrate will exhibit a decrease in value to dilution effect.
 Lipemic samples >+3 should be ultracentrifuged and the analysis performed on the supernatant.
 Do not use Sodium Fluoride because it inhibits urease activity, causing erroneous results.
 Refer to Literature References 3 and 4 for other interferences caused by drugs and diseases.

If the instrument fails calibration or controls do not measure within acceptable range when Diamond Diagnostics products are used, Diamond Diagnostics suggests the following:
 Verify that the reagents and internal calibrators used to standardize the instrument are correct for the instrument, have adequate expiration, and do not contain visually evident contamination.
 Follow the procedures delineated within the Operator’s Manual listed under Troubleshooting.
 Ensure that all appropriate Maintenance Procedures, as listed in the Operator’s Manual, have been performed.
 If problems still exist, contact Diamond Diagnostics’ Technical Service Department.

The names of various manufacturers, their instruments, and their products referred to herein may be protected by trademark or other law, and are used herein solely for purpose of reference. Diamond Diagnostics expressly disclaims any affiliation with them or sponsorship by them.

Useable Range:							
	<p>The Conductivity Rate Method used in the BUN assay on the Synchron CX® Systems provided the following useable range: Serum/ Plasma/ Urine (CX Delta only) – 3-150 mg/dL</p> <ul style="list-style-type: none"> - On the CX® Delta the Over Range Detection And Correction (ORDAC) function will accommodate samples with BUN concentrations between 150 to 300 mg/dL. - Samples with concentrations exceeding the high end of the analytic range should be diluted with saline and reanalyzed. 						
Functional sensitivity:							
	<p>Evaluated on dilutions of serum samples made from a starting serum of an approximate concentration of 13.7mg/dL; dilutions were 1:3, 1:5, 1:11 and a 0 mg/dL. Samples were tested as 4 samples per run over 5 calibrated runs.</p> <ul style="list-style-type: none"> - The lowest level where a CV of less than 20% was achieved was with the dilution at an expected value of 3 mg/dL BUN which recovered as: <table style="margin-left: 20px; border: none;"> <tr> <td style="padding-right: 20px;">3 mg/dL with Diamond reagent</td> </tr> <tr> <td>5 mg/dL with Beckman reagent</td> </tr> </table> 	3 mg/dL with Diamond reagent	5 mg/dL with Beckman reagent				
3 mg/dL with Diamond reagent							
5 mg/dL with Beckman reagent							
Reference Value:							
	<p>The values given as reference^{3,4,5,6} ranges are intended only as a guide rather than as absolute values. No adjustment has been made for age, sex, or dietary differences that may exist. Each lab should establish a reference range based on their patient population:</p> <table style="margin-left: 20px; border: none;"> <tr> <td style="padding-right: 20px;">Sample Type</td> <td>SI Units</td> </tr> <tr> <td>Serum or Plasma</td> <td>6-20 mg/dL</td> </tr> <tr> <td>Urine</td> <td>7-16 g/24hrs</td> </tr> </table>	Sample Type	SI Units	Serum or Plasma	6-20 mg/dL	Urine	7-16 g/24hrs
Sample Type	SI Units						
Serum or Plasma	6-20 mg/dL						
Urine	7-16 g/24hrs						

Performance Characteristics

Performance data was collected on CX®3 Delta and CX®9ALX Systems

Method Comparison:	<p>Diamond Bun Reagent was compared to Beckman Reagent using spiked or diluted serum samples run in triplicate and tested with each reagent, Diamond Bun Reagent and Beckman Bun Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained: $Diamond = 0.892 \times Beckman + 1.22 \text{ mg/dL}$ $r^2 = 0.997 \quad n=66 \quad \text{range} = 4 \text{ to } 87 \text{ mg/dL} \quad S_{(y,x)} = 1.30 \text{ mg/dL}$</p> <p>Diamond Bun Reagent was compared to Beckman Reagent using spiked or diluted urine controls run in triplicate and tested with each reagent, Diamond Bun Reagent and Beckman Bun Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained: $Diamond = 0.923 \times Beckman + 1.05 \text{ mg/dL}$ $r^2 = 0.996 \quad n=42 \quad \text{range} = 4 \text{ to } 100 \text{ mg/dL} \quad S_{(y,x)} = 1.70 \text{ mg/dL}$</p>
Recovery to Expected Values:	<p>Was evaluated for each matrix, serum and urine. Dilutions of the respective matrices were made and measured with Diamond and Beckman reagent. % Recovery = (Measured/expected) x100 was calculated for both Diamond and Beckman. Diamond and Beckman exhibited similar recoveries across the range of values in all matrices. See table 1 below.</p>
Precision:	<p>Commercially available control sera and urine controls were assayed using both Diamond and Beckman Creatinine Reagents on Synchron CX® Delta System. Samples were run for 20 days, 2 runs per day, 2 observations per run. See table 2 below.</p>

1

Matrix	Range of Conc. Expected, mg/dL	Reagent	Range of average % Recovery	Overall Mean Recovery
Serum	81.0 – 7.5 mg/dL	Diamond	82 – 106	97
		Beckman	94 – 111	105
Urine	100 – 10.5 mg/dL	Diamond	85.7 – 103.8	97.6
		Beckman	93.3 – 110	101.2

2

	N	Test Mean mg/dL	within run sd	Total sd
Serum Control 1	80	13	0.7	0.7
Serum Control 2	80	48	1.3	3.1
Urine Control 1	80	58	1.2	1.9
Urine Control 2	80	63	1.9	2.6

Literature References

1. Chin, W.T., and W. Kroontie, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757 – 1760 (1961)
2. Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
3. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
4. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
5. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Použitá literatura

- Chin, W.T., and W. Kroontie, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757 – 1760 (1961)
- Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
- Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
- Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
Pour plus de renseignements contactez:
Diamond Diagnostics - Support Technique +1-508-429-0450

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- Chin, W.T., and W. Kroontie, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757– 1760 (1961)
- Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
- Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
- Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Diamond Diagnostics, 333 Fiske Street, Holliston, MA.
For Technical Assistance call:
Diamond Diagnostics Technical Services at 1-508-429-0450

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Sample Type	SI Units
Serum or Plasma	6-20 mg/dL
Urine	7-16 g/24hrs

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Literatur Referenz

1. Chin, W.T., and W. Kroontje, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757 – 1760 (1961)
2. Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
3. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
4. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
5. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Irodalmi referenciák

- Chin, W.T., and W. Kroontje, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757 – 1760 (1961)
- Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
- Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
- Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Referências Literatura:
<ol style="list-style-type: none"> Chin, W.T., and W. Kroontie, "Conductivity Method for Determination of Urea", <i>Anal. Chem.</i> 33:1757 – 1760 (1961) Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", <i>Clin. Chem.</i> 17:644 (1971) Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994) Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001) SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Serviciul Technic Diamond Diagnostics 1-508-429-0450

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1				
Matrix	Range of Conc. Expected, mg/dL	Reagent	Range of average % Recovery	Overall Mean Recovery
Serum	81.0 – 7.5 mg/dL	Diamond	82 – 106	97
		Beckman	94 – 111	105
Urine	100 – 10.5 mg/dL	Diamond	85.7 – 103.8	97.6
		Beckman	93.3 – 110	101.2

2				
	N	Test Mean mg/dL	within run sd	Total sd
Serum Control 1	80	13	0.7	0.7
Serum Control 2	80	48	1.3	3.1
Urine Control 1	80	58	1.2	1.9
Urine Control 2	80	63	1.9	2.6

Referencias Literarias

1. Chin, W.T., and W. Kroontie, "Conductivity Method for Determination of Urea", *Anal. Chem.* 33:1757– 1760 (1961)
2. Paulson, G., R. Ray, and J. Sternberg, "A Rate-Sensing Approach to Urea Measurement", *Clin. Chem.* 17:644 (1971)
3. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
4. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
5. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA