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
Intended Use: Diamond Glucose Reagent is intended for the quantitative determination of glucose in serum, plasma, cerebrospinal fluid and urine on Beckman Synchron CX® & CX® Delta Systems.

Summary And Principle: This product is intended to serve as a functional equivalent to pre-existing material distributed by the Original Equipment Manufacturer (OEM).

Clinical Significance: Glucose measurements are used in the diagnosis and treatment of carbohydrate metabolism disorders including diabetes, mellitus, neonatal hypoglycemia, and pancreatic islet cell carcinoma.³

Methodology: Glucose concentration is determined by an oxygen rate method employing a Beckman Oxygen Electrode. Electronic circuits determine the rate of oxygen consumption, which is directly proportional to the concentration of glucose in the sample.¹

Chemical Reaction Scheme: When a 10 micro liters sample is injected into a reaction cup containing the glucose reagent, the glucose undergoes changes according to the following reaction:

Glucose from the sample reacts with the oxygen in the reagent in the presence of glucose oxidase (GO). The oxygen, which is consumed in direct proportion to the glucose, is monitored using the oxygen electrode, and the rate of depletion is a measure of the amount of glucose in the sample.

Reagents: Glucose Reagent, BK-443355D

Containing: Glucose Reagent – 500mL is an aqueous solution containing 590 U/mL Glucose Oxidase, 10% denatured alcohol, 0.04 mol/L Potassium Iodide, ammonium molybdate, buffer, and preservative.
Also required on the CX® for the determination of glucose: 1-10 liter bottle of Wash Solution (BK-443335D – reconstituted for 250mL concentrate); Synchron Calibration Standards; Appropriate Controls.

For in vitro diagnostic use only

Cautions: Exercise normal laboratory precautions. If contact occurs with skin, rinse affected area with water. If contact with eyes occurs, immediately rinse with copious amount of clean water or eye rinse. In cases of accidental ingestion, contact a physician immediately.

Storage and Stability: Glucose Reagent stored at 2°C to 8°C is stable until the expiration date indicated on each bottle. Glucose Reagent stored at room temperature (ie as when on the instrument) is stable for 30 days or until expiration date, which ever is sooner. Do no use if frozen.

Procedure

Reagent Preparation: Prior to use, allow the reagent to warm to ambient temperature. An 18 – 25°C water bath may be used to warm reagent.

Quality Control: Diamond Diagnostics suggests the use of commercially available control material with results assayed for the instrument used. Controls should be run at Normal and Abnormal levels. Diamond Diagnostics suggests measuring controls before patient samples are run and following instrument maintenance.

Limitations

Serum, Plasma or Urine are the sample of choice, whole blood samples are not recommended^{4,5}. Compatible anticoagulants or Chemical additives with this method^{5,6}:

Anticoagulants / Acceptable Level		Anticoagulants / Acceptable Level		Anticoagulants / Acceptable Level	
Sodium Citrate*	3.5 mg/mL	Ammonium Heparin	45 Units/mL	Sodium Heparin	45 Units/mL
EDTA	4.0 mg/mL	Lithium Heparin	45 Units/mL	Lithium Iodoacetate	1.5 mg/mL

*Specimens collected in a tube containing liquid Sodium Citrate will exhibit a decrease in value to dilution effect.
Lipemic samples >+3 should be ultracentrifuged and the analysis performed on the supernatant.
Refer to Literature References 2,3,4 and 5 for other interferences caused by drugs and diseases.

If the instrument fails calibration or controls do not measure within acceptable range when Diamond Diagnostics products are used, Diamond Diagnostics suggests the following:
Verify that the reagents and internal calibrators used to standardize the instrument are correct for the instrument, have adequate expiration, and do not contain visually evident contamination.
Follow the procedures delineated within the Operator's Manual listed under Troubleshooting.
Ensure that all appropriate Maintenance Procedures, as listed in the Operator's Manual, have been performed.
If problems still exist, contact Diamond Diagnostics' Technical Service Department.

Useable Range:	The Oxygen Rate Method used in the glucose assay on the Synchron CX® Systems provided the following useable range: Serum/ Plasma– 0.6 to 25 mmol/L or 10 to 450 mg/dL; Urine/ CSF 25 to 50 mmol/L or 450 to 900 mg/dL -On the CX® Systems the Over Range Detection And Correction (ORDAC) function will accommodate samples with Glucose concentrations between 450 to 900 mg/dL. -Samples with concentrations exceeding the high end of the analytic range should be diluted with saline and reanalyzed.												
Functional sensitivity:	Evaluated on dilutions of serum samples made from a starting serum of an approximate concentration of 42.5mg/dL; dilutions were 1:3, 1:5, 1:11 and a 0 mg/dL. Samples were tested as 4 samples per run over 5 calibrated runs. - The lowest level where a CV of less than 20% was achieved was with the dilution at an expected value of 9 mg/dL Glucose which recovered as: <table><tr><td>10 mg/dL with Diamond reagent</td></tr><tr><td>8 mg/dL with Beckman reagent</td></tr></table>	10 mg/dL with Diamond reagent	8 mg/dL with Beckman reagent										
10 mg/dL with Diamond reagent													
8 mg/dL with Beckman reagent													
Reference Value:	The values given as reference^{3,4,5} ranges are intended only as a guide rather than as absolute values. No adjustment has been made for age, sex, or dietary differences that may exist. Each lab should establish a reference range based on their patient population: <table><tr><td><u>Sample Type</u></td><td><u>Units mmol/L</u></td><td><u>SI Units</u></td></tr><tr><td>Serum or Plasma</td><td>3.9 - 5.8</td><td>70 - 105</td></tr><tr><td>Urine</td><td>0.06 - 0.8</td><td>1 - 15</td></tr><tr><td>CSF</td><td>2.2 – 3.9</td><td>40 - 70</td></tr></table>	<u>Sample Type</u>	<u>Units mmol/L</u>	<u>SI Units</u>	Serum or Plasma	3.9 - 5.8	70 - 105	Urine	0.06 - 0.8	1 - 15	CSF	2.2 – 3.9	40 - 70
<u>Sample Type</u>	<u>Units mmol/L</u>	<u>SI Units</u>											
Serum or Plasma	3.9 - 5.8	70 - 105											
Urine	0.06 - 0.8	1 - 15											
CSF	2.2 – 3.9	40 - 70											

Performance Characteristics	
Performance data was collected on CX®3 Delta and CX®9ALX Systems	
Method Comparison:	Diamond Glucose Reagent was compared to Beckman Reagent using spiked or diluted serum samples run in triplicate and tested with each reagent, Diamond Glucose Reagent and Beckman Glucose Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained: Diamond = 1.038 x Beckman – 2.31 mg/dL $r^2 = 0.998$ n=50 range = 0 to 900 mg/dL $S_{(y,x)} = 9.00$ mg/dL Diamond Glucose Reagent was compared to Beckman Reagent using spiked or diluted urine controls run in triplicate and tested with each reagent, Diamond Glucose Reagent and Beckman Glucose Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained: Diamond = 1.022 x Beckman + 1.067 mg/dL $r^2 = 0.998$ n=57 range = 1 to 293 mg/dL $S_{(y,x)} = 2.85$ mg/dL Diamond Glucose Reagent was compared to Beckman Reagent using spiked or diluted CSF run in triplicate and tested with each reagent, Diamond Glucose Reagent and Beckman Glucose Reagent. Recoveries of individual observations were compared by least squares regression and the following statistics were obtained: Diamond = 1.014 x Beckman – 0.920 mg/dL $r^2 = 0.997$ n=36 range = 8 to 118 mg/dL $S_{(y,x)} = 1.68$ mg/dL
Recovery to Expected Values:	Was evaluated for each matrix, serum, urine and CSF. Dilutions of the respective matrices were made and measured with Diamond and Beckman reagent. % Recovery = (Measured/expected) x100 was calculated for both Diamond and Beckman. Diamond and Beckman exhibited similar recoveries across the range of values in all matrices. See table 1 below.
Precision:	Commercially available control sera and urine controls were assayed using both Diamond and Beckman Creatinine Reagents on Synchron CX® Delta System. Samples were run for 20 days, 2 runs per day, 2 observations per run. See table 2 below.

1

Matrix	Range of Conc. Expected, mg/dL	Reagent	Range of average % Recovery	Overall Mean Recovery
Serum	835 – 30 mg/dL	Diamond	83.3 – 115.4	103
		Beckman	83.3 – 117.5	103
Urine	300 – 10 mg/dL	Diamond	88 – 120	101
		Beckman	87 – 110	97
CSF	100 - 25mg/dL	Diamond	88 – 98	94
		Beckman	88 – 100	95

2

	N	Test Mean mg/dL	within run sd	Total sd
CSF Control 1	80	59	0.8	3.8
CSF Control 2	80	30	1.5	3.0
Serum Control 1	80	89	2.1	6.0
Serum Control 2	80	308	2.2	21.7
Urine Control 1	80	45	1.0	2.3
Urine Control 2	80	287	2.4	14.4

Literature References

1. Kadish, A.H. and D.A. Hall, "A New Method for the Continuous Monitoring of Blood Glucose by Measurement of Dissolved Oxygen", *Clin Chem.* 11:869 (1965)
2. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
3. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
4. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA
5. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

Manufacturer and Product Information

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.

For Technical Assistance call:

Diamond Diagnostics Technical Services at 1-508-429-0450

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Literature References

- Kadish, A.H. and D.A. Hall, "A New Method for the Continuous Monitoring of Blood Glucose by Measurement of Dissolved Oxygen", *Clin Chem.* 11:869 (1965)
- Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
- Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA
- SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Diamond Diagnostics - Support Technique +1-508-429-0450

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Irodalmi referenciák

1. Kadish, A.H. and D.A. Hall, "A New Method for the Continuous Monitoring of Blood Glucose by Measurement of Dissolved Oxygen", *Clin Chem.* **11:869** (1965)
2. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
3. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
4. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA
5. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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6. Kadish, A.H. and D.A. Hall, "A New Method for the Continuous Monitoring of Blood Glucose by Measurement of Dissolved Oxygen", *Clin Chem.* 11:869 (1965)
7. Tietz, N. W., ed., "Clinical Guide to Laboratory Tests", 2nd Edition, W.B. Saunders Company, Philadelphia (1994)
8. Henry, John Bernard, "Clinical Diagnosis and Management by Laboratory Methods", 20th Edition, W.B. Saunders Company, Philadelphia (2001)
9. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA
10. SYNCHRON CX®3 Delta Clinical System Operating Instructions, Beckman Instruments Inc., Brea CA

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Serviciul Technic Diamond Diagnostics 1-508-429-0450

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