Manufacturer and Product Information	
Diamond Diagnostics, 333 Fiske Street, Holliston, MA. For Technical Assistance call: Diamond Diagnostics Technical Services at 1-508-429-0450	
Intended Use:	Used in conjunction with Calibrator 1 and 3 to calibrate Na, K, Cl, Ca, tCO ₂ , Glucose, BUN, & Creatinine on Beckman Synchron CX Delta, CX®CE, and CX9 ALX Systems.
Summary And Principle:	This product is intended to serve as a functional equivalent to pre-existing material distributed by the Original Equipment Manufacturer (OEM). Diamond Diagnostics calibrating standards have defined electrolyte concentrations that provide internal calibration points against which samples are measured.
Reagents:	Calibrator 2, BK-465908D, 6 x 25 mL
Containing:	Calibrator 2 is an aqueous solution containing: Na 180 mmol/L; K 10 mmol/L; Cl 180 mmol/L; Ca 14 mg/dL; Creatinine 8 mg/dL; Glucose 50 mg/dL; Urea Nitrogen 80 mg/dL; tCO ₂ 0 mmol/L; buffer; surfactant and preservative.
For <i>in vitro</i> diagnostic use only	
Values Assignment:	Na/K, & Cl values are referenced to an aqueous standard made with corresponding analyte NIST (National Institute of Standards and Technology) material (919a/918a, 919a respectively). Ca, Glucose, Creatinine, & Urea N values are referenced to a Calibrator (Beckman 465909 or a previous lot of Diamond Calibrator BK-465909) on a CX Delta or CX®3, except Ca.
Cautions:	Exercise normal laboratory precautions. If contact occurs with skin, rinse affected area with water. If contact with eyes occurs, immediately rinse with copious amount of clean water or eye rinse. In cases of accidental ingestion, contact a physician immediately.
Stability:	Product stability is listed on the product label. The product should not be used beyond this date. Store upright at 2° - 8 °C. After opening the calibration standard it is stable for 1 month if stored at Room Temperature (25°C) or until the expiration date, which ever is sooner. Do not Freeze. Benzoic acid crystals may form as the result of storage at lower temperatures. Crystals may be redissolved by warming or maintaining bottle at 18 – 26°C for 24 to 48 hrs.

Procedure	
Procedure:	The product is manufactured in a ready to use form. It is intended to serve as a direct replacement to pre-existing materials distributed by the OEM. For a detailed description of the use of this reagent, refer to the Instrument's Operator Manual.
Quality Control:	Diamond Diagnostics suggests the use of commercially available control material with results assayed for the instrument used. Controls should be run at Normal and Abnormal levels. Diamond Diagnostics suggests measuring controls before patient samples are run and following instrument maintenance.

Limitations	
Limitations:	If the instrument fails calibration or controls do not measure within acceptable range when Diamond Diagnostics products are used, Diamond Diagnostics suggests the following:
Verify that the internal calibrators used to standardize the instrument are correct for the instrument, have adequate expiration, and do not contain visually evident contamination.	
Follow the procedures delineated within the Operator's Manual listed under Troubleshooting.	
Ensure that all appropriate Maintenance Procedures, as listed in the Operator's Manual, have been performed.	
If problems still exist, contact Diamond Diagnostics' Technical Service Department.	

Calibrator 2

Beckman Synchron CX Delta, CX®CE & CX9 ALX Systems

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Diamond Diagnostics Technical Services at 1-508-429-0450

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