Manufacturer and Product Information

Diamond Diagnostics, 333 Fiske Street, Holliston, MA.

For Technical Assistance call:

Diamond Diagnostics Technical Services at 1-508-429-0450

Contents:
1 – pCO₂ Electrode
1 – Cleaning Solution (RD-943905E)
1 – Cuvette

Intended Use: Use with Radiometer ABL 5xx Series analyzers.

Storage: Store the pCO₂ electrode at room temperature (18-25°C) until you are ready to use it in your system.

Precautions: Following pCO₂ electrode replacement or membrane change, automatic calibrations are changed to the following Sequence:
1. First 4 hours, alternate in 30 minute intervals Gas Cal 1 and Gas Cal 2, ending in Gas Cal 2.
2. Next 20 hours, every hour Gas Cal 1, every second hour alternate with Gas Cal 2.

Procedure

Installation:

1. Unpack the new pCO₂ electrode. Preparation requires 25 hours.
2. Add 5 drops of Electrode cleaner, RD-943905E, into the Cuvette provided with the electrode.
3. Soak the electrode in the vial for 1 hour.
4. Remove the electrode and rinse with plenty of distilled or deionized water.
5. Empty the Cuvette and rinse with distilled or deionized water. Fill with distilled/deionized water up to one third the height of the Cuvette.
6. Soak the electrode in distilled/deionized water for 24 hours. The electrode may drift excessively if soaked for a shorter time.
7. Remove the cover of the temperature controlled electrode unit.
8. Press the pCO₂ amplifier housing down, pull the metal tab outward, and lift the amplifier housing.
9. Remove the old pCO₂ electrode.
10. Remove the newly prepared electrode from the Cuvette with water. Gently shake excess water from the electrode.
11. Obtain a box of pCO₂ membranes in jackets. Remove the foil cover.
12. Press the electrode firmly into the jacket.
13. Remove the jacketed electrode from the box.
14. Look for bubbles between the membrane and the electrode. If bubbles are observed, gently move the electrode up and down in the electrolyte to release the air bubbles. Press the electrode back into the jacket again.
15. Using a dry clean cloth or paper towel, thoroughly wipe and dry the gold contacts of the electrode.
16. Place the newly prepared electrode in the correct section of the measuring chamber.
17. Gently press the amplifier housing down onto the electrode.
18. Record in the maintenance schedule that the pCO₂ electrode has been replaced.
19. If all maintenance steps have been performed, place the cover over the temperature controlled electrode unit and press **Restart Analyzer** key below the display in scheduled maintenance menu

Warranty

Warranty = 12 months: Warranty is calculated from date of shipment from Holliston, MA.

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