



Mission Control™ Blood Gas and Electrolyte Control - Level 1

LOT R11E124-EU
2014/04

Expected Ranges Chart

Blood Gas/ISE Analyzer	pH			pCO ₂ mmHg			pO ₂ mmHg			Na ⁺ mmol/L			K ⁺ mmol/L			Ca ⁺⁺ mmol/L			Cl ⁻ mmol/L			Li ⁺ mmol/L		
	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max
Diamond PROLYTE									116	113	119	1.98	1.89	2.07				77	74	80				
Diamond SMARTLYTE, GEMLYTE									112	109	115	1.94	1.85	2.03	2.26	2.09	2.43	78	74	82	0.27	0.25	0.29	
Eschweiler Combiline	7.110	7.082	7.138	88	82	94	77	69	85	118	114	122	2.04	1.95	2.13	2.19	2.03	2.35	84	80	88	0.46	0.43	0.49
Eschweiler CombisysII	7.110	7.082	7.138	85	79	91	75	67	83	118	114	122	2.00	1.91	2.09	2.20	2.04	2.37	88	84	92	0.46	0.43	0.49
Eschweiler ECOLYTE	7.110	7.082	7.138							118	114	122	2.00	1.91	2.09	2.20	2.04	2.37	88	84	92	0.46	0.43	0.49
Eschweiler ECOSYSII	7.110	7.082	7.138	85	79	91	75	67	83															
Fresenius Ionometer									114	111	117	1.98	1.89	2.07	1.94	1.79	2.09							
Medica EasyLyte Na/K, Na/K/Cl, Na/K/Li, Na/K/Cl/Li, Na/K/pH/Ca	7.105	7.077	7.133						121	117	125	2.07	1.98	2.16	2.26	2.09	2.43	84	80	88	0.31	0.29	0.33	
Radiometer ABL 5	7.10	7.072	7.128	82	77	87	64	57	71															
Radiometer ABL 50, 500, 505, 510, 520	7.098	7.070	7.126	81	76	86	79	70	88	120	116	124	1.92	1.83	2.01	1.97	1.82	2.12						
Radiometer ABL 555	7.143	7.114	7.172	68	64	72	79	70	88	120	116	124	1.92	1.83	2.01	1.97	1.82	2.12						
Radiometer ABL 600, 610, 620, EML-100	7.106	7.078	7.134	85	79	91	70	62	78	120	116	124	1.92	1.83	2.01	1.97	1.82	2.12	79	75	83			
Radiometer ABL 700	7.106	7.078	7.134	85	79	91	68	61	75															
Roche/AVL 990, 995	7.086	7.058	7.114	86	80	92	73	65	81															
Roche/AVL 9110, 9140	7.106	7.078	7.134							112	109	115	2.11	2.02	2.20	2.11	1.95	2.27						
Roche AVL 9120, 9130										114	111	117	2.02	1.93	2.11				85	81	89			
Roche/AVL 9180, 9181										110	107	113	1.94	1.85	2.03	2.12	1.96	2.28	78	74	82	0.29	0.27	0.31
Roche/AVL Cobas b 121	7.136	7.107	7.165	82	77	87	55	49	61	121	117	125	2.12	2.02	2.22	1.88	1.74	2.02	87	83	91			
Roche/AVL Cobas b 221	7.126	7.097	7.155	82	77	87	28	25	31	121	117	125	2.12	2.02	2.22	1.88	1.74	2.02	87	83	91			
Roche/AVL Compact Series	7.086	7.058	7.114	86	80	92	73	65	81															
Siemens/Bayer 248	7.102	7.074	7.130	79	74	84	55	49	61															
Siemens/Bayer 348	7.113	7.085	7.141	80	75	85	67	60	74	116	113	119	2.20	2.10	2.30	1.85	1.71	1.99	88	84	92			
Siemens/Bayer 614, 634, 644, 654, 664	7.13	7.101	7.159							119	115	123	1.92	1.83	2.01	1.96	1.81	2.11	84	80	88	0.27	0.26	0.29
Siemens/Bayer 840, 845, 850, 855, 860, 865	7.127	7.098	7.156	78	73	83	66	59	73	111	108	114	1.99	1.90	2.08	1.93	1.79	2.07	81	77	85			
Siemens/Bayer RapidPoint 400, 405	7.126	7.097	7.155	88	82	94	60	53	67	116	113	119	1.82	1.74	1.90	1.82	1.68	1.96	80	76	84			

<p>for in vitro Diagnostic Use In vitro Diagnosticum Usage in vitro Использование in vitro Указан Аппарат для использования in vitro In vitro diagnosticum 器具 体外診断用 для использования в диагностике in vitro</p>	<p>European Conformity CE-Konformitätskennzeichnung Conformité aux normes européennes Conformidad europea Conformidade com as normas europeias Europäisch übereinstimmende 符合欧洲 Европейская Адекватность</p>	<p>Temperature Limit Temperaturlimit Limite de température Limite de temperatura Limite de temperatura Temperaturgrenze 温度限制 Температурные ограничения</p>	<p>Consult Instructions for Use Gebrauchsanweisung beachten Consulte la notice d'emploi Consulte las instrucciones de uso Consulte as instruções de utilização Beri brugervejledning 参考使用说明书 Рекомендации по применению</p>	<p>Lot Number Chargen-Nr. Número de lote Número de lote Número de lote Batch number 批号 Номер серии</p>	<p>Use by (YYYY-MM-DD) Verwendbar bis (JJJJ-MM-TT) Date de péremption (AAAA-MM-JJ) Usar hasta el (AAAA-MM-DD) Utilizar até (AAAA-MM-DD) Anvend før (AAAA-MM-DD) Использовать для (год-месяц-день рождения) 有效期至</p>	<p>Manufactured by Hergestellt von Fabriqué par Fabricado por Fremstillet af 製造者 оизводитель</p>	<p>Authorized Representative Bevollmächtigter Représentant agréé Representante autorizado Representante autorizado Autoriseret repræsentant 授權代表 Санкционированный представитель</p>	<p>Catalog Number Katalognummer Número de catalogue Número de catálogo Número de catálogo Katalognr. 产品编号 Номер каталога</p>
--	--	--	---	--	---	--	--	--