



Blood Gas/ISE Analyzer	pH			pCO <sub>2</sub> mmHg			pO <sub>2</sub> mmHg			Na <sup>+</sup> mmol/L			K <sup>+</sup> mmol/L			Ca <sup>++</sup> mmol/L			Cl <sup>-</sup> mmol/L			Li <sup>+</sup> mmol/L				
	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max	Mean	Min	Max		
Diamond PROLYTE												157	153	162	6.87	6.56	7.18				115	110	121	2.38	2.24	2.53
Diamond SMARTLYTE, GEMLYTE												152	148	157	6.42	6.13	6.71	0.51	0.44	0.58	115	110	120	2.05	1.93	2.17
Diamond SMARTLYTE PLUS												151	147	156	6.79	6.48	7.09	0.46	0.39	0.53	115	110	120	2.26	2.12	2.40
Diamond UNITY												154	149	158	6.84	6.53	7.15				120	115	126			
Erba Mannheim, EC 90												177	172	182	7.37	7.04	7.70	0.53	0.45	0.61	135	129	141			
Eschweiler Combiline	7.687	7.656	7.718	19.8	18.5	21.1	79	75	84	160	155	165	7.34	7.01	7.67	0.48	0.41	0.55	118	113	124	2.32	2.18	2.46		
Eschweiler Combisys II	7.687	7.656	7.718	19.8	18.5	21.1	79	75	84	160	155	165	7.34	7.01	7.67	0.48	0.41	0.55	121	116	127	2.32	2.18	2.46		
Eschweiler ECOLYTE												159	154	164	7.34	7.01	7.67	0.48	0.41	0.55	121	116	127	2.32	2.18	2.46
Eschweiler ECOSYS II	7.687	7.656	7.718	19.8	18.5	21.1	79	75	84																	
Horiba Yumizen E100												151	147	156	6.79	6.48	7.09	0.46	0.39	0.53	115	110	120	2.26	2.12	2.40
IL Iyte	7.701	7.670	7.732									164	159	168	7.30	6.97	7.63	0.83	0.71	0.95	123	118	129	2.37	2.23	2.52
Fresenius Ionometer												158	154	163	7.19	6.87	7.52	0.61	0.52	0.70						
Medica EasyLyte Na/K, Na/K/Cl, Na/K/Li, Na/K/Cl/Li, Na/K/pH/Ca	7.701	7.670	7.732									158	153	163	7.30	6.97	7.63	0.83	0.71	0.95	123	118	129	2.37	2.23	2.52
Radiometer ABL 5	7.64	7.61	7.67	23.3	21.8	24.8	80	76	84																	
Radiometer ABL 50, 500, 505, 510, 520, 555	7.706	7.675	7.737	18.8	17.5	20.0	85	81	90	168	163	173	7.50	7.16	7.84	0.68	0.58	0.78								
Radiometer ABL 600, 610, 620, EML-100	7.703	7.673	7.734	16.8	15.7	17.9	86	82	91	155	151	160	7.04	6.73	7.36	0.91	0.78	1.04	117	112	123					
Roche/AVL 990, 995	7.684	7.653	7.715	18.1	17.0	19.3	87	83	92																	
Roche/AVL 9110, 9140	7.684	7.653	7.714									153	149	158	6.36	6.08	6.65	0.66	0.56	0.76						
Roche AVL 9120, 9130												153	149	158	6.36	6.08	6.65				125	119	131			
Roche/AVL 9180, 9181												152	148	157	6.40	6.11	6.69	0.53	0.45	0.61	115	110	120	2.15	2.02	2.28
Roche/AVL Cobas b 121	7.674	7.643	7.705	17.1	16.0	18.2	76	72	80	167	162	172	7.30	6.98	7.63	0.63	0.54	0.73	126	120	132					
Roche/AVL Cobas b 221	7.663	7.633	7.694	18.1	16.9	19.3	76	72	80	165	160	170	7.30	6.98	7.63	0.63	0.54	0.73	124	118	129					
Roche/AVL Compact Series	7.688	7.658	7.719	17.1	16.0	18.2	85	81	90																	
Siemens/Bayer 248	7.623	7.592	7.653	24.0	22.4	25.5	76	72	81																	
Siemens/Bayer 348	7.652	7.621	7.682	20.7	19.4	22.1	87	82	92	165	160	170	6.45	6.16	6.74	0.57	0.48	0.65	116	111	121					
Siemens/Bayer 614, 634, 644, 654, 664	7.71	7.68	7.74							158	153	163	7.01	6.70	7.33	0.61	0.52	0.70	123	118	129	2.22	2.09	2.36		
Siemens/Bayer 840, 845, 850, 855, 860, 865	7.669	7.638	7.700	23.6	22.0	25.1	82	78	87	152	148	157	6.72	6.42	7.02	0.53	0.45	0.60	115	110	120					
Siemens/Bayer RapidPoint 400, 405	7.733	7.703	7.764	17.3	16.2	18.4	78	74	83	170	165	175	7.60	7.26	7.94	0.59	0.51	0.68	123	117	128					



For In Vitro Diagnostic Use  
In Vitro Diagnostikum  
Usage In Vitro  
Para Uso Diagnóstico In Vitro  
Usar apenas em diagnóstico in vitro  
Til in vitro diagnostisk  
仅供体外诊断使用



European Conformity  
CE-Konformitätskennzeichnung  
Conformité aux normes européennes  
Conformidade europea  
Conformidade com as normas europeias  
Europäischer Konformitätskennzeichen  
符合歐



Temperature Limit  
Temperaturlimit  
Limite de température  
Limite de temperatura  
Temperaturgrenze  
測和温度限制



Consult Instructions for Use  
Gebrauchsanweisung beachten  
Consulte le notice emploi  
Consulte las instrucciones de uso  
Consulte as instruções de utilização  
Benyttningsanvisning  
參閱說明書使用



Lot Number  
Charge-Nr.  
Número de lote  
Batchnummer  
批号



Use by (YYYY-MM-DD)  
Verwendbar bis (JJJJ-MM-TT)  
Date de péremption (AAAA-MM-JJ)  
Usar hasta el (AAAA-MM-DD)  
Utilizar até (AAAA-MM-DD)  
Anvend til (AAAA-MM-DD)  
效期(YYYY-MM-DD)



Manufactured by  
Hergestellt von  
Fabriqué par  
Fabricado por  
Fabricado por  
Fremstillet af  
\*\*\*製造



Authorized Representative  
Bevollmächtigter  
Représentant agréé  
Representante autorizado  
Autorisierter repräsentant  
授權的代理



Catalog Number  
Katalognummer  
Número de catálogo  
Número de catálogo  
Katalog  
产品编号