



# DIAMOND®

DIAGNOSTICS



# CLINICHECK™

定值血清基质冻干生化控制品



## CliniCheck™ 定值生化质控品

### 详细说明

基于人体血清  
 保质期: 两年(2-8°C)  
 冷冻复溶的稳定性: 20天(大多数分析物)  
 开启后稳定性: 7天(2-8°C下, 大多数分析物)  
 不需要特殊的稀释液, 各种酶和CO2的冷冻复溶稳定性良好  
 定值于超过20个仪器

### ANALYTES

Acetaminophen	Cholesterol (LDL)	LDH
Acid Phosphatase (Total)	Cholesterol (Total)	Lipase
Albumin	Cholinesterase	Lithium
Alkaline Phosphatase (ALP)	CO2	Magnesium
Alpha-1-Antitrypsin (1)	Complement C3	Phosphorus
Alpha-Fetoprotein (1)	Complement C4	Potassium
Alpha-Hydroxybutyrate Dehydrogenase (aHBDH)	Creatine Kinase (CK)	Protein (Total)
ALT/SGPT	Creatinine	Salicylate
Amylase	Digoxin	Sodium
Apolipoprotein A-1	GGT	T3 (Free)
AST/SGOT	Glucose	T3 (Total)
Bilirubin (Direct)	Haptoglobin	T4 (Free)
Bilirubin (Total)	IgA	T4 (Total)
Calcium	IgG	Transferrin
Carcinoembryonic Antigen (1)	IgM	Triglycerides
Ceruloplasmin (1)	Iron	TSH
Chloride	Iron (TIBC)	Urea
Cholesterol (HDL)	Iron (UIBC)	Urea Nitrogen
	Lactate (Lactic Acid)	Uric Acid

### 包装

DD-93001 水平1	12 x 5 mL
DD-93002 水平2	12 x 5 mL
DD-93012 双水平	6 x 2 x 5mL
DD-93100 校正液和稀释液	6 x 2 x 3mL