


DIAMOND®

DIAGNOSTICS



MISSION®

适用于血气, ISE, 血氧以及代谢物分析系统的质量控制品



Mission CONTROL™

水基质质控品

- 用于测试 pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, tCO₂
- 室温下保质期: 36个月
- 开启后60秒内性状稳定

适用仪器

Abbott I-Stat
Idexx
Instrumentation Laboratory
ITC (formerly Diametrics)
Medica
Nova Biomedical
OPTI Medical
Radiometer
Roche Diagnostics (formerly AVL)
Siemens Diagnostics (formerly Bayer)

Mission TRINITY™

含有染料和各种代谢物的水基质质控品

- 定值于 pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, Glucose, Lactate, tHb, sO₂, MetHb, COHb, HHb, O₂Hb
- 冷藏下保质期: 36个月
- 开启后60秒内性状稳定

适用仪器

Instrumentation Laboratory
Radiometer
Roche Diagnostics (formerly AVL)
Siemens Diagnostics (formerly Bayer)

Mission COMPLETE™

水基质CVM线性质控品


- 定值于 pH, pCO₂, pO₂, Na⁺, K⁺, Ca⁺⁺, Cl⁻, Li⁺, Glucose, Lactate
- 冷藏下保质期: 36个月
- 开启后60秒内性状稳定
- 含有5个不同颜色标识的安瓶, 可覆盖仪器重要的测量范围

适用仪器

Abbott I-Stat
Instrumentation Laboratory
ITC (formerly Diametrics)
Nova Biomedical
OPTI Medical
Radiometer
Roche Diagnostics (formerly AVL)
Siemens Diagnostics (formerly Bayer)

Mission WebSTAT™

- 安全, 简便地使用基于网络的统计数据软件
- 从月到期, 从年到期, 同年范围内质控参考
- 每月提供质控 Levy-Jennings 曲线图
- 完整客户还可在线使用线性回归图



The names of manufacturers, their instruments, and their products referred to herein may be protected by trademark or other law and are used herein solely for purpose of reference. Diamond Diagnostics expressly disclaims any affiliation with them or sponsorship by them. For current regulatory status on products within this catalog, please contact your sales representative.