

DIAMOND[®]

DIAGNOSTICS



CareLyte[®] PLUS

The Choice for Veterinary Laboratories



CareLyte® PLUS

The Choice for Veterinary Laboratories

Veterinary Options

Feline, Canine, Bovine, Equine, Swine, Ovine or Open

Sample

Whole Blood, Serum, Plasma, Urine or Dialysate

Sample Size

95 µl

Number of Parameters

4 (Simultaneously)

Detection Range

Blood Na⁺: 40-200 mmol/L
Dialysate K⁺: 1.7-15 mmol/L
Cl⁻: 50-200 mmol/L
Ca⁺⁺: 0.3-5.0 mmol/L
Li⁺: 0.2-5.5 mmol/L

Urine Na⁺: 3-300 mmol/L
K⁺: 5-120* mmol/L
Cl⁻: 15-300 mmol/L

*60-120 requires additional dilution

Reproducibility

Typical Within Run (n=30)

Blood, Serum, Plasma, Dialysate

Na⁺: CV ≤ 1% (120-160mmol/L)
K⁺: CV ≤ 2% (2.5-6 mmol/L)
Cl⁻: CV ≤ 2% (85-130 mmol/L)
Ca⁺⁺: SD ≤ 0.02 mmol/L (0.8-1.5 mmol/L)
Li⁺: SD ≤ 0.02 mmol/L (0.4-1.3 mmol/L)

(Lithium not typically measured in dialysate samples)

Urine

Na⁺: CV ≤ 5% (100-250 mmol/L)
K⁺: CV ≤ 5% (10-60 mmol/L)
Cl⁻: CV ≤ 5% (100-250 mmol/L)

(Calcium and Lithium are not typically measured in urine samples)

Display Resolution

Na⁺: 1 mmol/L or 0.1 mmol/L
K⁺: 0.1 mmol/L or 0.01 mmol/L
Cl⁻: 1 mmol/L or 0.1 mmol/L
Ca⁺⁺: 0.01 mmol/L or 0.001 mmol/L
Li⁺: 0.01 mmol/L or 0.001 mmol/L

This analyzer is produced and distributed by Diamond Diagnostics Inc. There is no affiliation, connection, sponsorship or association between Diamond Diagnostics Inc. and Roche Diagnostics or Idexx Laboratories.



Analysis Time

45 sec without printout (72 per hour)
60 sec with printout (60 per hour)

Data Storage

1000 Patient Results
500 Level QC Results (3 levels)

Calibration

2 Point every 4 hours
On Demand
1 Point after every Measurement

User Input

Touch Screen
USB Keyboard / Barcode Scanner (optional)

Output

Color 5" display 800 x 480 pixels
Graphic 16 column printer
LAN / RS232 LIS Communication

Features

Wifi, Bluetooth, 4 USB ports, LAN

Ambient Conditions

15-32°C, Room temperature (60-90°F)
<85% humidity

Languages

English, Chinese, French, German, Italian, Korean
Japanese, Portuguese, Russian, Spanish, Polish

Power

100-240V ~ 50/60 Hz (self adjusting)
1.6 A max, 50 Watts max

Size & Weight

13.2 x 12.4 x 12" or 335 x 315 x 295 mm
14 lbs or 6 kgs

DIAMOND®
D I A G N O S T I C S

333 Fiske Street, Holliston, MA 01746
Tel: 508-429-0450 | Fax: 508-429-0452
www.diamonddiagnostics.com

SOP23-0101F REV04 08/11/17