



SmartLyte[®] *Plus*

Alat Analisa Elektrolit Tercanggih

SmartLyte[®] Plus



Parameter	Jarak	Reproduksibilitas ¹	Resolusi
Na ⁺	40-200 mmol/L	CV ≤ 1% (120-160mmol/L)	0.1 mmol/L
K ⁺	1.7-15 mmol/L	CV ≤ 1.5% (2.5-6 mmol/L)	0.01 mmol/L
Cl ⁻	50-200 mmol/L	CV ≤ 1% (85-130 mmol/L)	0.1 mmol/L
Ca ⁺⁺	0.3-5.0 mmol/L	SD ≤ 0.02 mmol/L (0.8-1.5 mmol/L)	0.001 mmol/L
Li ⁺	0.2-5.5 mmol/L	SD ≤ 0.03 mmol/L (0.4-1.3 mmol/L)	0.001 mmol/L
Urin ² Na ⁺	3-300 mmol/L	CV ≤ 5% (100-250 mmol/L)	0.1 mmol/L
K ⁺	5-120 mmol/L ³	CV ≤ 5% (10-60 mmol/L)	0.01 mmol/L
Cl ⁻	15-300 mmol/L	CV ≤ 5% (100-250 mmol/L)	0.1 mmol/L

1-Tipe dalam pemeriksaan (n=30) darah, serum, plasma 2-Kalsium dan litium bukan tipe ukuran dan sampel urin
3-(60-120) memerlukan tambahan pengencer

Daya

100-240V ~ 50/60 Hz (bisa disesuaikan)
1.6 A maksimum, 50 Watts maksimum

Ukuran Dan Berat

13.2 x 12.4 x 12 in or 335 x 315 x 295 mm
14 lbs or 6 kgs

Kondisi lingkungan

Suhu ruang: 15-32°C/60-90°F
Kelembaban <85%

Tidak dijual di Amerika Serikat