



DIAMOND[®]

DIAGNOSTICS



CLINICHECK[™]

Lyophilized Serum based Assayed Chemistry Control



CliniCheck™ Assayed Chemistry Control

SPECIFICATIONS

Human serum based
 2 year shelf life at 2–8°C
 20 day frozen reconstituted stability for most analytes
 7 day open-vial stability at 2–8°C for most analytes
 Requires no special diluent and has good frozen reconstituted stability for enzymes and CO₂.
 Assayed for over 20 instruments.

ANALYTES

Acetaminophen	Cholesterol (LDL)	LDH
Acid Phosphatase (Total)	Cholesterol (Total)	Lipase
Albumin	Cholinesterase	Lithium
Alkaline Phosphatase (ALP)	CO ₂	Magnesium
Alpha-1-Antitrypsin (1)	Complement C3	Phosphorus
Alpha-Fetoprotein (1)	Complement C4	Potassium
Alpha-Hydroxybutyrate Dehydrogenase (aHBDH)	Creatine Kinase (CK)	Protein (Total)
ALT/SGPT	Creatinine	Salicylate
Amylase	Digoxin	Sodium
Apolipoprotein A-1	GGT	T3 (Free)
AST/SGOT	Glucose	T3 (Total)
Bilirubin (Direct)	Haptoglobin	T4 (Free)
Bilirubin (Total)	IgA	T4 (Total)
Calcium	IgG	Transferrin
Carcinoembryonic Antigen (1)	IgM	Triglycerides
Ceruloplasmin (1)	Iron	TSH
Chloride	Iron (TIBC)	Urea
Cholesterol (HDL)	Iron (UIBC)	Urea Nitrogen
	Lactate (Lactic Acid)	Uric Acid

PACKAGING

DD-93001 Level 1	12 x 5 mL
DD-93002 Level 2	12 x 5 mL
DD-93012 Dual Level	6 x 2 x 5mL
DD-93100 Calibrator and Diluent	6 x 2 x 3mL